

# T Cell Health Annexin V Assay (Non-Lyse)


## INSTRUCTIONS FOR ANALYSIS

### INTENDED USE


These instructions are intended for analyzing data obtained with the Accellix T Cell Health Annexin V Assay (Non-Lyse) run on the Accellix Platform. The data can be analyzed with FlowJo™ Software (BD Biosciences) using the **gating template Accellix T Cell Health Annexin V Assay (NL) Gating Template.wsp** provided with the assay or accessible in the Customer Resources section of [accellix.com](http://accellix.com). For additional FlowJo Software instructions, refer to [Accellix Data Analysis Guide for FlowJo Software](#), also available in the Customer Resources section on [accellix.com](http://accellix.com). Equivalent software can also be used for data analysis.

### ANALYSIS PROTOCOL


1. Retrieve the **events.fcs** files for analysis in FlowJo Software or an equivalent analysis software.
2. Download the **gating template** and use it to pull each FMO/FS pair into the workspace for analysis.
3. See the recommended gating strategy in **figure 1**.
  - a. Use this gating strategy to identify the early-stage apoptosis cells (Q3: Annexin-V+, Viability-) and necrotic cells (Q2: Annexin-V+, Viability+).
  - b. Identify T cell subsets, including viable CD45+ leukocytes (viable CD45+), CD3+ T cells, CD4+ T helper cells, and CD8+ cytotoxic T cells.
4. Example FlowJo plots for the T Cell Health Annexin V (NL) assay using the recommended gating strategy are provided in **figure 2**. Use a conservative cells (tight) gate to identify the cells in the CD45/FSC plot. Be sure to exclude debris (low forward scatter events).
5. Apply the FMO for the Viability/Annexin V plot to each FS in the pair.
  - a. Use a 5% contour setting to place the vertical and horizontal gates for the Viability/Annexin V plots.
  - b. Place the horizontal gate at the density minimum in the FMO (**figure 3**).
  - c. Apply the FMO (full stain minus one) gates to the FS (full stain) run (**figure 3**).
6. Place the gates at the density minimum for the Viability/FSC plot (**figure 4**).



**Figure 1.** Proposed gating strategy for Accellix T Cell Health Annexin V (Non-Lyse) Assay.




**Figure 2.** Cryopreserved PBMCs thawed and stained with the T Cell Health Annexin (V) NL assay



Example 1

Example 2

**Figure 3.** Examples of FMO and FS gate placement.



**Figure 4.** Examples of viability plot gate placement.